

# BED SETTEE MULTY

## Designer

Claude Brisson

## Technical characteristics

Multy offers an intermediate position half-way between settee and bed: a chaise longue with ample leg support. Frame in multi-ply and panels of particleboard, with steel mechanism and birch slats (thickness 9 mm). Seat in either polyether foam (28kg/m<sup>3</sup> - 3.0kPa) or 9 cm thick high resilience polyurethane Bultex foam, density 38 kg/ m<sup>3</sup>. Back in polyether foam (28kg/m<sup>3</sup> - 3.0kPa). Back in polyether foam (28kg/m<sup>3</sup> - 3.0kPa). Lumbar cushions are supplied as follows: 2 x W60 for small bed-settee without arms, and 2 x W78 for large bed-settee without arms; 1 x W78 for small bed-settee with arms, and 2 x W60 for large bed-settee with arms. Please take these details into consideration when ordering additional/subsequent lumbar. The armrests (polyether foam 28 kg/m<sup>3</sup> - 3.0 kPa + polyester quilting 200 g/m<sup>2</sup>) or bedside tables (MDF top diameter 43 cm - thickness 2 cm) pivot to allow the bed-settee to open. MULTY is easy to operate, and offers excellent comfort both as a bed (with polyether mattress size : 153x194x8, 119x194x8) and as a settee. When unfolded, its depth is 200 cm for the bed-settee W 155. Mattress has an ecru-coloured cover in 50% polyester/50% cotton (dry cleanable), which has been treated with "HealthGuard" (antibacterial & fungicidal). A foam sheet is supplied to place between the slats and the seat



## DIMENSIONS

when folded. For instructions on use, please see the card which comes with the bed. MULTY has been specially designed for use as a bed on an everyday basis. Optional extras: lumbar cushions and arms. Covers of lumbar cushions and armrests are removable.

**More information on**

[www.ligne-roset.com](http://www.ligne-roset.com)

© Ligne Roset 2023